

Microplate Seals

Film and Foil Microplate Heat Seals

A comprehensive range of film and foil microplate heat seals, ideal for use with the Vitl VTS and Vitl Micro^{TS}



Description

Vitl stock a comprehensive range of film and foil microplate heat seals, ideal for (but not limited to) use with the Vitl VTS Microplate Heat Sealer and Vitl Micro^{TS} Microplate Heat Sealer.

All Vitl microplate seals are supplied in packs of 100 sheets, each measuring 125mm x 80mm.

Reels of seals are available on request.

Technical Specifications

	Peelable?	Pierceable?	Temperature range	Seals Polypropylene?	Seals Polyethylene?	Seals Polystyrene?	Seals COC?	Color code	Catalogue Number
20µm Heat Sealing Foil	X	✓	-20 - +120°C	✓	✓	✓	X	Green	V901002
38µm Heat Sealing Foil	✓	✓	-20 - +110°C	✓	✓	✓	X	Blue	V901003
70µm Heat Sealing Foil	✓	X	-80 - +80°C	✓	✓	X	X	Red	V901001
75µm Clear Sealing Film	X	✓	-80 - +80°C	✓	✓	✓	✓	75µm	V901005
85µm Heat Sealing Foil	✓	X	-200 - +110°C	✓	X	X	X	Black	V901004
105µm Clear Sealing Film	✓	X	-80 - +180°C	✓	✓	✓	✓	105µm	V901006

Related Products



VTS
Variable
Temperature
Sealer
V902001



Micro^{TS}
Microplate
Heat Sealer
V903001

Vitl 
LIFE SCIENCE SOLUTIONS
A member of the ITL group

T: +86(021)58585463

vitlproducts.com